

WOW

UNE MARQUE DE BIOTECH DENTAL

# SCAN·MATE

BY ARMOR DENTAL

YOUR INTRAORAL  
SCANNING  
PARTNER

**Comfortable retraction, control, and efficiency  
that produces quality images.**

- ✓ Improving digital work flow for intraoral scanning.
- ✓ Error-free scans from the beginning.
- ✓ Latex-free and fully autoclavable instrument up to 60 cycles.



GET YOUR SCAN RIGHT THE FIRST TIME.  
YOU'RE WELCOME!



## PROFESSIONAL ISOLATION, RETRACTION INSTRUMENT FOR INTRAORAL SCANNING AND ANY RETRACTION NEEDS:

- Ultimate retraction instrument for all dental procedures or exams.
- Non-reflective device for enhanced scanning for single restoration to full arches.
- Ideal for all intraoral scanning needs – crowns, bridges, implants, and orthodontics.
- Comfortable retraction for scanners, lasers, intraoral cameras, and air polishers.
- Durable ergonomic handle to improve control and quality of image.
- Soft, bendable tip provides smooth retraction and easier access for hard to reach areas.

---

## SCAN•MATE IS CHANGING EVERYTHING FOR DIGITAL WORK FLOW

- Reduces scanning time for dental team and results in cleaner scans.
- Improves patient experience and first impression chairside.
- More efficient, results in high quality scans.
- No more tissue issues when scanning.
- Helps the dental team improve their scanning skills, freeing up dentists' time to do other procedures.

---

## COMFORTABLE FOR CLINICIAN AND PATIENT

- Scan-Mate can be used for scanning and all retraction needs.
- No more fingers or mirrors or other isolation products needed.
- Patented bendable neck and ergonomic handle make it easy to scan images with comfort – bite image, 2nd posterior molar, lower distal, unruly tongue, tight lips.
- Reduces clean-up time after scanning.
- Latex-free.

---

## TO AUTOCLAVE SCAN-MATE

- Steam autoclave at 250F to 273°F (121°C to 134°C) for at least 20 minutes as per ADA guidelines. (Time may vary depending on size of load and autoclave.)
- Autoclavable 60 cycles and then replace after 1 month.
- Unique bendable handle and soft tip improves the patient experience.
- Scan•Mate's adjustable neck will eventually weaken after a month.
- Do not keep over the recommended autoclave limit to ensure product effectiveness.



BIOTECH DENTAL

Manufacturer: Armor Dental Corporation - Agent: MedPass SAS - Distributor: Biotech Dental. Class I device, CE.

Manufacturer WOW: Condor SAS - Distributor: Biotech Dental. Class I medical device for digital dentistry. CE.

All brands are under the responsibility of their respective manufacturers. Read the instructions carefully. No refund by social security.

Images are for representation purpose only. Do not throw in public areas.

Biotech Dental - S.A.S with a capital of 28 866 417 € - RCS Salon-de-Provence: 795 001 304 - SIRET: 795 001 304 00018 - VAT N °: FR 31 79 500 13 04.